

Rubber Fab Technologies Group's  
**COMPONENT  
PRODUCTS**

Pumps, Valves, Fillers  
and System Equipment



## Our Commitment Continues

Rubber Fab Technologies Group has **successfully met our customer's goals** with **quality products** that meet the **highest sanitary standards and purity criteria**. Our sanitary gaskets have provided **superior surface finish and elastomer purity** for use in production facilities of fine pharmaceutical, biotechnology, high-purity water, injectables, food and beverage products. In fact, Rubber Fab has made **many significant advancements in gasket materials and manufacturing** over the past few years and has recently applied this knowledge to the production of our **Component Parts for pumps, valves, fillers and system equipment**.


Rubber Fab has made the **same commitment to purity and excellence** in the manufacturing of our **Component Parts** as we have our line of gaskets. Like our other products, our component parts are **color coded for material identification** and come with **lot and batch numbers that provided traceability** at all times.

If it's **valve stems and balls, centrifugal pump parts and kits or filler parts**, you can be assured of our continuous **commitment for a superior product**. Rubber Fab Technologies Group **will meet your requirements and your expectations**.

### Valve Stems and Balls:

- New and exchange Valve Stems
- Check Valve Replacement Balls
- PTFE enhanced
- Superior surface finishes
- Available in buna, EPDM and viton
- FDA and 3A compliant
- All stem varieties are available: compression, divert and shut-off

## Rubber Fab Replacement Parts for Series Centrifugal Pumps



### Centrifugal Pump Kits for:

- 114, 216 & 218/328 pumps
- Available in buna, EPDM and Viton elastomers

### Centrifugal Pump Parts:

- Backplates
- Carbon seals
- Casing gaskets & O-rings
- Silicon carbide seals
- Cups
- Springs
- Stub shafts
- Impeller pins

### Filler Parts:

- Rolling diaphragms
- Vacuum cups
- Suction cups
- U cups
- Screens
- Perforated plates
- O-Rings
- Gaskets
- Couplings



Casing Gaskets, Carbons, O-Rings

## Metal Detectable Elastomers

Another major **material advancement** from Rubber Fab Technologies Group is in the area of **metal detectable elastomers**.

A common problem with component and filler parts is wear and tear. When rubber breaks off it can migrate through your system and into your product. **Searching for and locating fragmented rubber parts is a costly, time consuming and inconclusive process** requiring expensive x-ray equipment, manual observation and an extensive maintenance program. Not locating a worn rubber piece can have an even costlier outcome.

Rubber Fab provides an **easier and cost effective solution with metal detectable elastomers**. By manufacturing any standard elastomer with a **3A, FDA approved** metal filler, a worn off piece of elastomer material **can now be located by an in-line metal detector**. This alarm enables your system to **instantly be shut down for the removal of contaminated product** and the **immediate replacement of the worn part**. Detectables **eliminate product recall, lowers product loss and decreases the risk of elastomers consumed in finished product**.

Detectables are available in buna, EPDM, FKM and silicone materials and are available in sanitary and cam lock gaskets, valve stems, check ball, o-rings and almost any molded product.



Metal detectable elastomers, if worn or broken off, can now be easily identified with in-line metal detection.

## Waukesha & VNE Check Balls

PART NUMBER	VALVE SIZE	ACTUAL O.D.	MATERIAL
RF 122-2-1.5	1.5"	2-1/8"	BUNA
RF 122-2E-1.5	1.5"	2-1/8"	EPDM
RF 122-2V-1.5	1.5"	2-1/8"	VITON
RF 122-2N-1.5	1.5"	2-1/8"	NYLON
RF 122-3-2	2"	2-5/8"	BUNA
RF 122-3E-2	2"	2-5/8"	EPDM
RF 122-3V-2	2"	2-5/8"	VITON
RF 122-3N-2	2"	2-5/8"	NYLON
RF 122-47-2.5	2.5"	3-5/16"	BUNA
RF 122-47E-2.5	2.5"	3-5/16"	EPDM
RF 122-47V-2.5	2.5"	3-5/16"	VITON
RF 122-47N-2.5	2.5"	3-5/16"	NYLON
RF 122-26-3	3"	4-1/4"	BUNA
RF 122-26E-3	3"	4-1/4"	EPDM
RF 122-26V-3	3"	4-1/4"	VITON
RF 122-26N-3	3"	4-1/4"	NYLON

## Tri-Clover Check Balls

PART NUMBER	VALVE SIZE	ACTUAL O.D.	MATERIAL
RF 45BY1.5-U	1.5"	1.615	BUNA
RF 45BY1.5-E	1.5"	1.615	EPDM
RF 45BY1.5-SFY	1.5"	1.615	VITON
RF 45BY1.5-NY	1.5"	1.615	NYLON
RF 45BY2-U	2"	2.115	BUNA
RF 45BY2-E	2"	2.115	EPDM
RF 45BY2-SFY	2"	2.115	VITON
RF 45BY2-NY	2"	2.115	NYLON
RF 45BY2.5-U	2.5"	3.490	BUNA
RF 45BY2.5-E	2.5"	3.490	EPDM
RF 45BY2.5-SFY	2.5"	3.490	VITON
RF 45BY2.5-NY	2.5"	3.490	NYLON
RF 45BY3-U	3"	4.115	BUNA
RF 45BY3-E	3"	4.115	EPDM
RF 45BY3-SFY	3"	4.115	VITON
RF 45BY3-NY	3"	4.115	NYLON



Rubber Fab Technologies Group • Limecrest Road  
 P.O. Box 626 • Andover, NJ 07821  
 phone: 973-579-2959 fax: 973-579-7275

[www.rubberfab.com](http://www.rubberfab.com)

## Valve Stems

RF19-571	
RF19-635	
RF19-1024	
RF19-561-10	
RF19-627	
RF19-1258	
RF19-1254	
RF19-1300	
RF19-1336	
RF19-627B-4	
RF19-627C-4	
RF19-1254-4	
RF19-1068B	
RF19-1308	
RF19-643	
RF19-647	
RF19-748	
RF19-1257	
RF19-1308 & 1308A	
RF19-1259	
RF19-821	
RF19-535R	
RF19-561-21	
RF19-572	
RF19-1003	
RF19-1276-1.5	
RF19-1276 (2, 2.5, 3)	
RF19-733	
RF19-1306	
RF19-1255A-1.5	
RF19-1255A-2	
RF19-1255A-2.5	
RF19-1255B-3	

Rubber Fab Technologies Group is a member of:

