



Rubber Fab Technologies Group

# GAUGE GUARD ISOLATOR GASKETS

Protect Critical Instrumentation  
With No Performance Loss



The Rubber Fab Gauge Guard Isolator gaskets **protect expensive diaphragm and instruments from damage without affecting instrument performance.** By combining a quality Rubber Fab gasket with a membrane, the Gauge Guard Isolator gasket **isolates hot or caustic saline solution from instrumentation.** A necessary and **cost effective** step in **assuring long life and accurate results from your instruments and gauges.** As some manufacturers will suggest the use of a gauge guard gasket, others note that the **failure to use a gauge guard gasket can void most instrument and gauge warranties.**




## All the Rubber Fab Benefits in a Gauge Guard Gasket

- Gasket membrane **protects instrumentation from caustic saline solutions**
- **Will not interfere with gauge operation or accuracy**
- Works with most industry standard instruments
- Help **extend the life** of gauges
- Available in:
  - FKM** - .010 membrane thickness
  - Platinum Silicone** - .010 membrane thickness
  - EPDM** - .010 membrane thickness
  - TUF-STEEL** - .005 membrane thickness
- Samples available upon request

## All Rubber Fab Gaskets Meet the Stringent Standards For Purity

- U.S. Pharmacopeia Class VI Certification
- Cytotoxicity Criteria
- Title 21 CFR 177.2600
- Title 21 CFR 177.1550
- USDA and 3-A Sanitary Standards
- Current Good Manufacturing Practices (CGMP)
- ASME-BPE Standards
- Animal Derived Ingredient Free




### FKM

1/2"	42GGI-SFY-050
3/4"	42GGI-SFY-075
1"	40GGI-SFY-100
1-1/2"	40GGI-SFY-150
2"	40GGI-SFY-200




### Platinum Silicone

1/2"	42GGI-PX-050
3/4"	42GGI-PX-075
1"	40GGI-PX-100
1-1/2"	40GGI-PX-150
2"	40GGI-PX-200



### TUF-STEEL

1/2"	42GGI-TS-050
3/4"	42GGI-TS-075
1"	40GGI-TS-100
1-1/2"	40GGI-TS-150
2"	40GGI-TS-200



### EPDM

1/2"	42GGI-E-050
3/4"	42GGI-E-075
1"	40GGI-E-100
1-1/2"	40GGI-E-150
2"	40GGI-E-200

## New From Rubber Fab Technologies Group



**Tuf-Flex**  
The World's Finest Gaskets



**V²B Gaskets**  
Tank Vent and Breaker System Gaskets



Rubber Fab Technologies Group  
Limecrest Road

P.O. Box 626 • Andover, NJ 07821

phone: 973-579-2959 fax: 973-579-7275

[www.rubberfab.com](http://www.rubberfab.com)

Rubber Fab Technologies Group is a member of:



International Association of  
Food Industry Suppliers



Distributed by: